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## Product Description

Visalys® CemCore is a dual-curing composite containing fluoride for adhesive cementation of indirect restorations and root posts. Moreover, Visalys® CemCore can be used to fabricate radiopaque core build-ups and build-up fillings.

Due to the Active-Connect-Technology (ACT) Visalys® CemCore includes an additional initiator system for the optimization of the adhesive bond.

Visalys® CemCore in the self-mixing 1:1 automix syringe is free of bisphenol A and its derivatives. BPA or BPA precursors are not used in the manufacturing process of this device.

Visalys® CemCore is easy to use and has high adhesive strength and good flowability with stable consistency.

Visalys® CemCore shall only be used in combination with self-etching Visalys® Tooth Primer. Visalys® Restorative Primer is recommended as the primer for the restoration surfaces.

Visalys® CemCore is offered in the shades Universal (A2/A3), Dark (A4), Opaque, Bleach and Translucent.

## Indications for Use

Permanent cementation of:

- crowns, bridges, inlays, onlays, veneers and adhesive bridges made of metal, ceramic, composite, hybrid and oxide ceramics
- root posts and indirect core build-ups

Composite for:

- adhesive core build-ups

## Patient Target Group

Patients being treated in the context of a dental procedure.

## Contraindications

Known allergies to methacrylic acid ester, benzotriazole compounds, benzoyl peroxide, camphorquinone, 4-methoxyphenol, or other constituents of Visalys® CemCore or Visalys® Tooth Primer (see Visalys® Tooth Primer instructions for use).

If there are no adequate means of ensuring drainage.

## Warnings

Allergic reactions to the product cannot be ruled out in sensitive persons. Sensitization through contact with skin/ mucous tissue is possible.

Materials containing eugenol can inhibit setting of composites. The use of hydrogen peroxide solutions can impair the adhesive strength. Retraction paste containing iron(III)sulfate can impair curing and adhesion of the material.

In the case of vital teeth, areas near the pulp may have to be treated with suitable pulp/dentin protection.

Keep out of reach for children.

## Cautions

Do not use after the expiration date.

Do not swallow material. If the material is inadvertently swallowed: Consult a medical doctor if problems arise or persist.

Take the appropriate protective measures (e.g. wear gloves or safety glasses). Avoid contact of skin and mucous tissue with non-cured material.

Avoid contact with the eyes. If there is accidental contact with the eyes: Rinse immediately and thoroughly with an eye wash or water. Consult a medical doctor if problems arise or persist.

Avoid contact with clothing, as permanent damage to clothes cannot be excluded.

Safety data sheets are available at [www.kettenbach.com](http://www.kettenbach.com).

## General Notes

Given refrigerator storage, the product has to have reached room temperature prior to use.

In case of contact with Visalys® Tooth Primer and direct incident radiation of ambient light, curing of Visalys® CemCore accelerates.

For layers thicker than the specified polymerization depth, it is necessary to wait until the self-curing process is complete in all cases.

Working times apply for products, which are processed at a temperature of 23°C/74°F and a relative humidity of 50%. Higher temperatures reduce the working time and lower temperatures extend these times.

## Application Instructions for Cementation

After removing the temporary appliance, carefully remove any residues of temporary cement from the core/cavity surface.

## Checking the Fitting Accuracy and Occlusion of the Restoration

Thin and brittle restorations can break under excessive pressure during try-in. For this reason, do not apply high pressure or only check the occlusion after permanent cementation.

With the available Visalys® CemCore Try In Pastes, the shade effect of the selected composite shade can be checked on the slightly moist core prior to final cementation. Apply the paste in the selected shade on the inside surfaces of the restoration and check the shade effect on the restoration. The shade of the Visalys® CemCore Try In Paste matches the shade of the cured cementation composite.

Remove the Visalys® CemCore Try In Paste from the restoration and core surface with water.

After the try-in, the surfaces of the restoration that may be contaminated with saliva, blood or silicone residues, for example, have to be cleaned with water spray or in the ultrasonic bath with distilled water and then dried with water/oil-free air. All the surfaces must be free of contamination.

In the case of restorations made of non-precious metal or oxide ceramics, there may be no cleaning with phosphoric acid.

A rubber dam is recommended for a dry operating field.

### **Pretreatment of the Restoration, the Root Post or the Indirect Core Build-Up**

Unless prescribed otherwise in the instructions for use for the restoration material, root post or material used for the indirect core build-up, the following procedure should be selected for pretreatment:

- 1.1 Adhesive surfaces of noble metal, metal and oxide ceramics (zirconium oxide and aluminum oxide ceramics)
  - Sand blasting of the surfaces with aluminum oxide powder ( $\leq 50 \mu\text{m}$ ). Observe the restoration material manufacturer's blasting parameters.
  - Then rinse off thoroughly with water spray or clean in the ultrasonic bath.  
DO NOT CLEAN WITH PHOSPHORIC ACID!
  - Dry with water / oil-free air.
- 1.2 Adhesive surfaces of silicate ceramic (feldspar ceramic, glass ceramics) and hybrid ceramics
  - Etch the surfaces with hydrofluoric acid according to the specifications in the instructions for use of the restoration material.
  - Then rinse off thoroughly with water spray or clean in the ultrasonic bath.
  - Dry with water / oil-free air.
- 1.3 Adhesive surfaces of composite
  - Sand blasting of the surfaces with aluminum oxide powder ( $\leq 50 \mu\text{m}$ ). Observe the restoration material manufacturer's blasting parameters.

- Then clean thoroughly with alcohol.
  - Dry with water / oil-free air.
2. Drop Visalys® Restorative Primer into a preparation tray. Light protection is not necessary. Close primer bottle carefully after each use.
  3. Apply Visalys® Restorative Primer with a disposable applicator or brush onto the adhesive surfaces and allow it to work in for 60 sec. Then blow out and dry with a gentle jet of water / oil-free air. Visalys® Restorative Primer is non-light curing.

The surfaces of some types of posts can be attacked by sandblasting or hydrofluoric acid. The specifications of the root post manufacturer on surface treatment must be observed.

### **Pretreatment of the Core / Cavity Surface (Enamel, Dentin, Metal, Composite)**

Visalys® CemCore may only be used together with Visalys® Tooth Primer. Visalys® Tooth Primer is a self-etching single-component primer and can therefore be used without additional phosphoric acid etching. Only with cementation of veneers and adhesive bridges and with uncut enamel is selective phosphoric acid etching necessary. Depending on the clinical conditions, in other cases too, selective phosphoric

acid etching may be optionally performed. Use the etching gel according to the instructions for use from the relevant manufacturer.

1. Following the try-in of the restoration, saliva, blood and any residues of the Visalys® CemCore Try In Paste must be removed from the adhesive surfaces of the dental hard substance and restorations surfaces using water spray and then keep the surfaces dry. The adhesive surfaces must be free of contamination. A rubber dam is recommended for a dry operating field.
2. Drop Visalys® Tooth Primer into a preparation tray. Light protection is not necessary. After applying the Visalys® Tooth Primer, start with the further processing immediately. Close primer bottle carefully after each use.
3. Apply Visalys® Tooth Primer to the entire adhesive surface of the tooth surface using a suitable applicator and rub in for 20 sec. Any excess can be removed with a cotton pellet and blown off with air. Then dry with a gentle jet of water / oil-free air. The prepared primer should be processed within 5 min., as constituents can evaporate. The Visalys® Tooth Primer is non light-curing.

For cementation on an implant abutment or on an indirect core build-up, the specifications of the respective material manufacturer must always be considered. Pretreat the entire

adhesive surfaces of the abutment/core build-up dependent on the material, as described under “Pretreatment of the restoration” in the Visalys® CemCore instructions for use, unless prescribed otherwise in the relevant material manufacturer’s instructions for use.

The emergence profiles and screw channels, as applicable, must be protected by suitable measures (e.g. wax, silicone covers).

### Application of the Cementation Composite

1. Apply Visalys® CemCore on the inside surfaces of the restorations and on the core or cavity surfaces, as applicable. In the case of application in the mouth, the processing time is reduced and the restorations must be placed within 2 min.
2. Place the restoration on the core or in the cavity.
3. Removing excess
- 3.1. **Method 1:** Remove after initial curing through brief exposure to light or self-curing until the “gel phase”

The excess is either polymerized with a commercially available polymerization lamp (1200 mW/cm<sup>2</sup>) for 2–3 sec. oral and vestibular as well as mesial and distal – i.e. at four points – or wait for the gel phase in self-cur-

ing. The excess, which is still soft-elastic, can be easily removed with a scaler.

#### 3.2. **Method 2:** Removing without initial curing

Remove the excess directly, e.g. with a cotton pellet, brush or foam pellet.

While removing the excess, the restoration has to be held in its final position. All excess, especially in difficult-to-reach approximal areas, has to be removed prior to final curing.

#### 4. Final curing

After removal of excess, the cement joint can be covered with a glycerin gel or polyethylene glycol paste, e.g. Visalys® CemCore Try In Paste, to avoid the formation of an inhibition layer. This can then be removed with water.

Then, polymerize each surface / cement joint for 10 seconds with light.

For non-translucent restorations, the final curing must be allowed to complete beneath the restoration (see the table for technical data).

Translucent restorations can also be cured in several layers according to the specified polymerization times dependent on the size of the light aperture.

## 5. Finishing

Rough cement joints can be finished and polished with the appropriate instruments.

### **Special Procedure for Cementation of Root Posts and Indirect Core Build-Ups**

1. Prepare the root canal, root post or indirect core build-up according to the instructions for use from the relevant manufacturer. If the manufacturer of the root post or indirect core build-up recommends silanization or pre-treatment with a primer for metal or oxide surfaces or a primer-adhesive / universal adhesive-activator combination, the Visalys® Restorative Primer can be used for this purpose.
2. Whetten the cleaned and dried surfaces of the root canal, the coronal dentin and the enamel parts, as applicable, completely with Visalys® Tooth Primer and rub in for 20 sec. Remove excess from the root canal with a paper tip and dry all surfaces with a gentle jet of water/oil-free air.
3. Apply Visalys® CemCore into the canal lumina and also on the root post. The end of the tip used should always remain in the material while applying. Do not use Lentulo spirals.

4. Immediately afterwards, insert the root post (the working time for Visalys® CemCore in the root canal is approx. 2 min.). It is recommended that any excess material is specifically distributed on the occlusal core surface.
5. To fix the post, the material can be light-cured with a suitable polymerization lamp for 10 sec.
6. Perform the build-up of the core as described under "Application instructions for core build-up".

### **Application Instructions for Core Build-Up**

1. Clean and dry the cavity.
2. Areas close to the pulp must be protected with the appropriate means.
3. Ensure a clean and dry operating field (rubber dam) during the entire application.
4. Selective phosphoric acid etching of the enamel may be undertaken optionally.
5. Drop Visalys® Tooth Primer into a preparation tray. Light protection of the primer is not necessary.
6. Apply Visalys® Tooth Primer to the entire adhesive surface of the core or the cavity using a suitable applicator and rub in for 20 sec. Any excess can be removed with a cotton pellet and blown off with air. Then dry with a jet of water /oil-free air.

- The prepared primer should be processed within 5 min., as constituents can evaporate.
- Place a matrix and apply Visalys® CemCore directly into the cavity, as applicable. The use of a transparent matrix makes it easier for light to enter in the case of optional light curing. Alternatively, in the case of very badly damaged teeth, a core build-up mold can be used. The mold is filled with Visalys® CemCore and then placed on the residual tooth substance. Visalys® CemCore can still be modeled within the processing time.
  - Wait for final curing (see Technical Data). Optionally, the material can be light-cured for 20 sec. to accelerate curing down to a depth of approx. 2.5 mm. Light-curing should be performed with a commercially available polymerization lamp (1200 mW/cm<sup>2</sup>). In case of a lower light intensity, prolong the light curing duration accordingly. Place the lamp as close to the material as possible.  
  
Deeper layers that are not accessible to light will cure in a purely chemical manner after approx. 5 min.
  - For layers thicker than the specified polymerization depth, it is necessary to wait until the self-curing process is complete in all cases.
  - After full curing, remove if applicable any matrices or core build-up formers.
  - Performing core preparation.

- Isolation of the core prior to fabrication of a temporary restoration.

In the case of self-curing, lower temperatures, e.g. of the product or during work in the anterior region, lead to prolongation of the curing time.

## Application

The paste is mixed and delivered directly from 9 g/5 ml-auto-mix syringe and the appropriate mixing tips. Only use original mixing tips (ML Ø 2.5 mm, blue) and Intraoral/Endo tips!

The following points have to be observed:

- Remove the syringe cap and prior to first attaching a mixing tip on the syringe, expel a small amount until both components are evenly delivered. Then discard the paste expelled.
- Attaching the mixing tip. Turn the cap 90° clockwise to lock in place. If applicable, mount an Intraoral/Endo tip locking on the mixer.
- The material requires a certain amount of time to flow out after the mixing tip has been attached. The material flow stops as soon as pressure is reduced on the syringe piston. Any material already cured in the mixing tip must not be extruded with force, because otherwise the syringe and the mixing tip may be damaged.

4. Dispose of used mixing tip, replace the syringe cap and check the outlet at the start of the next application.

## Cleaning and Disinfection

The risk of contamination of the automix syringe can be reduced by covering the syringe with a disposable protective sleeve.

Clean the syringe with a disposable wipe and standard media.

Disinfection of the syringe in accordance with national guidelines.

Observe standard practice hygiene measures.

## Storage Conditions

Store dry, protected from sunlight. The product has to be kept in the refrigerator at 2°C/36°F–8°C/46°F before first use.

## Shelf-Life After First Use

- a) For continued storage in the refrigerator:  
Until the end of the specified expiration date
- b) For continued storage at room temperature:  
Use within six months

Avoid extreme temperatures. Given refrigerator storage, the product has to have reached room temperature prior to use.

Do not use after the expiration date.

## Constituents

Main components of Visalys® CemCore are approx. 42 vol. % inorganic fillers in the size range 0.2–20 µm including ytterbium fluoride. The polymer base consists of aliphatic dimethacrylates.

## Classification

ISO 4049, type 2, class 3.

## Disposal

The empty automix syringes/mixing tips/applicator tips and cured material can be disposed in the household rubbish.

## Trademarks

Visalys® CemCore, Visalys® CemCore Try In Paste, Visalys® Tooth Primer and Visalys® Restorative Primer are registered trademarks of Kettenbach GmbH & Co. KG.

## Note

Not all Kettenbach products are available in all countries.

For dental professional use only. Caution: Federal (U.S.A.) Law restricts this device to sale by or on the order of a dentist.

If serious incidents occur that are related to the product, please immediately contact the manufacturer and the competent authority.

The Summary of Safety and Clinical Performance can be found in the European database on medical devices (EUDAMED – <https://ec.europa.eu/tools/eudamed>).

## Limitation of Liability

Except where prohibited by law, Kettenbach GmbH & Co. KG will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the legal basis asserted, including warranty, contract, negligence, or strict liability.

The information provided for Kettenbach products is based on comprehensive research and experience in application technology. Results are furnished to the best of our knowledge, subject to technical changes within the framework of product development. However, users must comply with and consider all recommendations and information in connection with any use.

Updated: November 21, 2024

Color stability meets the requirements of ISO 4049.

\* in a light wavelength range 400–500 nm

\*\* Aluminum has a radiopacity equivalent to dentin. Therefore, 1 mm of a material that has a radiopacity equivalent to 1 mm of aluminum has a radiopacity equivalent to 1 mm of dentin, and 2 mm of aluminum is equivalent to the tooth enamel.

## Technical Data

Working time (23°C/74°F)	approx. 4 minutes
Working time (intraoral)	approx. 2 minutes
Exposure time for initial curing (tack curing)	2–3 seconds per polymerization point
Time for purely chemical curing including working time (intraoral)	approx. 5 minutes
Light curing* (cementation)	10 seconds per surface/cement joint
Light curing* (core build-up)	20 seconds
Radiopacity**	approx. 2.5 mm Al
Depth of cure Translucent	approx. 2.5 mm
Depth of cure Universal (A2/A3)	approx. 2.0 mm
Depth of cure Bleach, Dark (A4)	approx. 1.5 mm
Depth of cure, Opaque	approx. 0.5 mm
Light intensity	1200 mW/cm <sup>2</sup>

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